

# Final Fill and Cryopreservation

Partnering with our customers to design complete, closed systems specific to their needs




**Challenge**

Filling final product into small aliquots can be time consuming and require open-transfer steps.

**Solution**

Disposables set, including overmolded assembly connected to small cryopreservation bags, allows for final filling to occur in one, closed step.

Concept



## Components

- VueLife® bags
- Overmolded C-Flex® with smooth inner bore
- Custom fluid transfer assembly
- Custom small volume KryoSure® cryopreservation bags

## Design

